# Vessel –on-a-chip: advanced in vitro model idis for cardiovascular research

Laura Vázquez-Vázquez<sup>1</sup>, Svitlana Palii<sup>2</sup>, Antón L. Martínez<sup>2,3</sup>, Lois Rivas-Meizoso<sup>1</sup>, Bastián Carnero<sup>1,4</sup>, Elias Ferreiro-Vila<sup>1</sup>, Sylvana V. Varela-Ballesta<sup>1</sup>, Nuria Chantada<sup>2</sup>, M. Teresa Flores-Arias<sup>4,</sup> José Manuel Brea<sup>2,3</sup>, Mabel Loza<sup>2,3</sup>, Bruno K. Rodiño-Janeiro<sup>1</sup>, Ezequiel Álvarez<sup>2,5,6.</sup>

1. BFlow S.L., Edificio Emprendia, Campus Vida s/n, 15706, Santiago de Compostela, Spain.

2. Department of Pharmacology, Pharmacy and pharmaceutical technology, Universidade de Santiago de Compostela (USC), 15782, Santiago de Compostela, Spain.

3. Drug Screening Platform/Biofarma Research Group, CIMUS Research Center, Universidade de Santiago de Compostela (USC), 15782, Santiago de Compostela, Spain.

4. Photonics4Life Research Group, Applied Physics Department, Faculty of Physics and Institute of Materials (iMATUS), Universidade de Santiago de Compostela, 15782 Santiago de Compostela, Spain.

5. Instituto de Investigación Sanitaria of Santiago (IDIS), Fundación IDIS, SERGAS, 15706, Santiago de Compostela.

6. CIBERCV, Madrid, Spain.


#### INTRODUCTION

Atherosclerosis is a chronic inflammatory disease of the arterial vessels that is behind of the main vascular pathologies as myocardial infarction and stroke. Atherosclerosis is a complex multifactorial process, but the physical factors such as flow-induced shear stress have an important role both in the onset and the progression of endothelial dysfunction. Traditional 2D endothelial models were useful to understand the biological factors that lead to atherosclerosis, but new in vitro models are needed to study the flow effects into the atherosclerosis mechanisms and the underlying therapeutic targets. Bflow models include relevant features such as circular section and a wide range of channels sizes In this work we present three models which are designed to study the following features:

- Vessel-on-a-chip with a bifurcation to characterize the local flow effects on the endothelial layer, including a model to functionalize with extracellular matrix (ECM).
- Atherosclerosis-on-a-chip with a bifurcation with an atherosclerotic plaque in one branch.
- Microvascular & Barrier chip with a region to stablish an extracellular matrix to perform assays of microvessels, angiogenesis and/or endothelial barrier.

## **METHODS**


Vessel-on-a-chip devices from BFlow were cultured with human umbilical vein endothelial cells. Different flow conditions were stablished with an external peristaltic pump. Different stainings were performed into the chips including live staining with AM-calcein, hoescht 33342, phalloidin and VE-cadherin immunostaining. Numerical simulations were performed to evaluate the precise flow effects into the chip.



A) Image of the chip before the experimental setup. B) Epifluorescence microscopy image of HUVEC stained with AM-calcein and hoescht 33342 covering the walls of the chip. C) Numerical simulation of the bifurcation showing the differential shear stress. D) 3D bifurcation reconstruction using confocal microscopy of HUVEC stained with phalloidin and VE-cadherin. E) HUVEC immunostained with VE-Cadherin antibody and secondary antibody Alexa 488.


On the left, image of a ECM Vessel-on-a-chip before the experimental setup. On the right, image of the loaded chip where the ECM is confined (blue) and the cell culture medium is filling the bifurcation (pink).

## **Microvascular & Barrier chip**



On the left, image of the extracelular matrix loading into the chip. On the right, schematic image of the sagital section of the chip indicating where the extracelullar matrix is confined and the complete

#### Atherosclerosis chip



On the left, image of the chip before the experimental setup. On the right, wall shear stress numerical

simulation of the atherosclerosis model, with one healthy branch and the flow alterations by the circular section of the lateral channels. atherosclerosis plaque in the other branch.

#### CONCLUSIONS

Bifurcated Vessel-on-a-chip reproduces the local flow effects into big arteries, being possible to functionalize these regions with ECM. Atherosclerosis chip reproduce the strong flow effects induced by the plaque. Microvascular & Barrier chip is a relevant model to study both microvascular formation and endothelial barrier function. In this work we presented four new models which will allow to improve our knowledge in endothelial biology, the cardiovascular pathophysiology and to design new therapies.

#### **BIBLIOGRAPHY**



Check the related papers here

# **ACKNOWLEDGEMENTS**

Laura Vazquez and Lois Rivas-Meizoso are supported by "INVESTIGO Program" contract from Xunta de Galicia. Sivtlana Palii is supported by Erasmus+ KA107 programme. Laura Vázquez-Vázquez, Lois Rivas-Meizoso, Bastián Bastián Carnero is supported by "Industrial PhD" contract no. 11 IN606D 2021 2604925 by GAIN/Xunta de Galicia. Elias Ferreiro-Vila is supported by "CONECTACOVID" Carnero, Elias Ferreiro-Vila, Sylvana V. Varela-Ballesta project DAVID by GAIN/Xunta de Galicia. Sylvana V. Varela-Ballesta is supported by "CONECTACOVID" project NANOFLUIDITEC by GAIN/Xunta de Galicia. Bruno K. Rodiño and Bruno K, Rodiño-Janeiro are employees of BFlow Janeiro is supported by "Torres Quevedo 2020" contract by Agencia Estatal de Investigación, Ministerio de Economía y Competitividad from Spain. The Project is S.L.. M. Teresa Flores-Arias and Ezequiel Álvarez are supported by "NEOTEC 2022" no. VOAC-ARPY SNEO-20222232 by Centro para el Desarrollo Tecnológico Industrial (CDTI), Spain, and "Agencia Estatal de Investigación, founders and partners of BFlow S.L. Ministerio de Economía y Competitividad" (Spain) under contract RTI2018-097063-B-100, by the Consellería de Educación Xunta de Galicia/FEDER (ED431B 2020/29 UE) and "la Caixa" Banking Foundation CaixaImpulse Validate 2021 Call, gran number CI20-00200.

## **CONFLICT OF INTEREST**

CARDIOCHUS